Effect of Cross-Linking on Mechanical Function in the Degenerate Nucleus Pulposus

An Nguyen

Wade Johannessen, Dawn M. Elliott

McKay Orthopaedic Research Laboratory University of Pennsylvania

SUNFEST Final Presentation, July 29, 2005



McKay Orthopaedic Research Laboratory

Motivation

- Low back pain
 - #2 reason Americans see their doctor
 - #2 reason for missing a workday
 - #1 cause of disability
- Intervertebral disc (IVD) degeneration




Intervertebral Disc

- Nucleus Pulposus (NP)
 - High proteoglycan content
 - Hydrated gel
- Annulus Fibrosus (AF)
 - High collagen content
 - Organized fibers




Intervertebral Disc Degeneration

Complex process



Proteoglycan content

- Decreases in NP with degeneration
- Changes in mechanical properties
 - Significant decrease in swelling pressure




Proteoglycan (Alberts, *Molecular Biology of the Cell*)



Cross-Linking

- Provide mechanical loading support
- Formed by reactions between amino side groups
- Cross-linking agent: Genipin











Intramolecuar

Study Hypothesis

Hypothesis: Cross-linking will restore mechanical function in the degenerate nucleus pulposus




ChABC Dose Finding Results





Genipin Cross-Linking Results





Next Step: Quantifying Cross-Linking

Spectrophotometry

- Measures amount of light absorbed by a substance
- Beer's Law: [Absorbance] = k[Concentration]
- Standards of known concentrations are used to make a standard curve






Ninhydrin Assay

- Spectrophotometry protocol for free amino acid group determination
- More cross-linking
 - Lower concentration of free amino acid groups
 - Lower absorbance measured



Preliminary Ninhydrin Data



Future Work

- Finalize ChABC dosage
- Test ChABC and genipin together
- Quantify genipin cross-linking

By understanding the effects of cross-linking, an alternative treatment for disc degeneration and low back pain may be developed.



Thank You



MCKAY ORTHOPAEDIC RESEARCH LAB

UNIVERSITY OF PENNSYLVANIA





McKay Orthopaedic Research Laboratory